


CLEANING A DRIP STAND


Drips stands are found to be frequently contaminated by pathogens which are able to survive for weeks or months causing healthcare associated infections. This shared piece of equipment is regularly moved between patients so is at a high risk of spreading germs if not decontaminated properly. Before cleaning please inspect for surface damage.

1/6




Perform hand hygiene. Wear recommended PPE.

2/6



Remove one wipe from the Universal pack.

4/6




Change wipe if it is dry or soiled and discard. Allow surfaces cleaned with Universal wipes to air dry.

3/6

1. TOP SURFACE


2. FRAME

3. WHEELS



Clean in order of: top surface, the frame and lastly wheels.

5/6



Remove all PPE and discard as per local policy. Perform hand hygiene.

6/6



As per local policy, add cleaned date and sign the indicator tape/note.