clinell®

Alcohol Wipes

Surface disinfection of non-invasive medical devices.



Rapid, high quality alcohol wipes

Ultra-quick drying, pharmaceutical grade 70% isopropyl alcohol wipes.

Disinfection wipes, which leave no smear, for surface disinfection of non-invasive medical devices.

Proven efficacy



High-quality durable wipes containing pharmaceutical grade 70% isopropyl alcohol. Our wipes are proven to have bactericidal and virucidal action.

Unique packaging



Wipes are dispensed one at a time, so the top of the next wipe is not dry, which commonly occurs with tub dispensers. This eliminates the infection control risk of dry contaminated wipes.

Less plastic. Less waste



Superior quality



Our flow wrap packaging is environmentally friendly as there is much less packaging to dispose of compared to standard alcohol tubs.

Clinell Alcohol Wipes contain 25gsm non-woven spun bond, which is wetter & thicker than competing alcohol wipes allowing for larger area disinfection with fewer wipes, saving time and money.

For the surface disinfection of non-invasive medical devices:

- Ultrasound probes
- Pulse oximeter probes
- Doppler probes
- Stethoscopes
- ECG machines

- Blood pressure machines
- Blood pressure cuffs
- Radiology equipment
- Incubators
- Monitors

Antimicrobial efficacy



Proven to kill at least 99.999% of microorganisms.

Including MRSA, *E. coli, Pseudomonas, Enterococcus, Candida* and *Aspergillus* according to EN1276, EN13727, EN14561, EN13624 & EN14562.

	Microorganism example	EN test
Bacteria	Acinetobacter baumannii	EN13727
	Enterococcus hirae	EN13727
		EN 16615
	Escherichia coli K12	EN13727
	Methicillin resistant Staphylococcus aureus	EN13727
	Pseudomonas aeruginosa	EN13727
		EN16615
	Staphylococcus typhimurium	EN13727
	Staphylococcus aureus	EN13727
		EN16615
Fungi	Aspergillus brasiliensis	EN13624
Mycobacteria	Mycobacterium avium (2 year aged)	EN14348
	Mycobacterium terrae (2 year aged)	EN14348
Viruses	Norovirus	EN14476
	Murine coronavirus	EN14476
	Vaccinia virus	EN14476
Yeast	Candida albicans	EN13624
		EN16615

Ordering information

Alcohol Wipes - Pack of	150		Transport unit view
ALCOHOL WIPES ALCHOL WIPES ALCHI MINISTER CHARLES ALCHI MINISTER A	Product code	CAW150L	1284.0
	Wipes per pack	150	
	NHSSC code (UK only)	VJT839	
	Expiry date	2 years	
	Unit per carton	18	
	Carton size	390 x 330 x 190 mm	
	Carton per layer	9	
	Layer height	1284 mm	
	Cartons per pallet	54	
	Classification	Class IIa Medical Device CE Mark	

For more information, speak to your GAMA Healthcare Area Manager or visit **www.clinell.com**



Use disinfectants safely. Always read the label and product information before use. Always follow the medical device manufacturer's cleaning and disinfection guidelines prior to use with Clinell Alcohol Wipes.

JBN220359

