## Clinell Peracetic Acid Wipes\* have been approved for use by the following manufacturers:

BRAND	EQUIPMENT
Abecca	Positioning Supports & Mattresses
Active Key GmbH & Co. KG	All Active Key MedicalKey™ products, including: AK-C4112, AK-C4412, AK-C7012 & AK-C7412
AGFA	DX-D products
Agile Medical	All products
BD (Bard)	Site~Rite® 8 Ultrasound System and Probe
Carrflex	Mattress Covers, Upholstery Tops & Associated Equipment
Carr Plastics	Techmaflex
Dartex Coating	Mattresses (with polyurethane face)
Fresenius Kabi - Agilia range	VP MC volumetric infusion pump SP MC syringe infusion pump SP TIVA TCI syringe infusion pump SP PCA syringe infusion pump
Fujifilm Healthcare	Ultrasound Probes and Accessories (Suitable for most, please see cleaning instruction QAD-Z3510 for full details)
Fujifilm Sonosite	Ultrasound Equipment (Edge, NanoMaxx, MicroMaxx, TITAN) iViz devices and C60v, L25v, L38v & P21v probe
Howard Wright	Mattresses

\*Previously Clinell Sporicidal Wipes.

Approved Manufacturers. © GAMA Healthcare Ltd. Version 14.

www.gamahealthcare.com/compatibility

JBN230355

**GAMA Healthcare Ltd.**, The Maylands Building, Maylands Avenue, Hemel Hempstead, Hertfordshire, HP2 7TG, UK. +44 (0)20 7993 0030 | info@gamahealthcare.com | www.gamahealthcare.com



Clinell Peracetic Acid Wipes\* have been approved for use by the following manufacturers:

BRAND	EQUIPMENT
ΗP	HP EliteBook 840 G6 Healthcare Edition HP EliteOne 800 G5 Healthcare Edition All-in-One HP Healthcare Edition HC270cr Clinical Review Display HP Healthcare Edition HC241 Clinical Review Monitor HP Healthcare Edition HC241p Clinical Review Monitor HP Healthcare Edition HC271 Clinical Review Monitor HP Healthcare Edition HC271p Clinical Review Monitor
Mindray	N Series Monitor
Modsel	All product ranges
OES Medical	Pressure Monitor and UAM Anaesthetic Machine
Orfit Industries	All related types of AIO 3.0 Low Density Head Supports Medium Density Blue Head Support
Panaz	Cadet Furniture
Parity Medical	Mobile Computer Carts
Philips	Ultrasound Equipment Intellivue Monitor (remaining residue must be wiped off) Ultrasound Transducers (Please email info@gamahealthcare.com for full list of compatible transducers)

\*Previously Clinell Sporicidal Wipes.

Approved Manufacturers. © GAMA Healthcare Ltd. Version 14.

www.gamahealthcare.com/compatibility

JBN230355

**GAMA Healthcare Ltd.**, The Maylands Building, Maylands Avenue, Hemel Hempstead, Hertfordshire, HP2 7TG, UK. +44 (0)20 7993 0030 | info@gamahealthcare.com | www.gamahealthcare.com



## Clinell Peracetic Acid Wipes\* have been approved for use by the following manufacturers:

BRAND	EQUIPMENT
Recticel	Mattresses
Sidhil	IQ Bed Frame
Siemens	Siemens RapidPoint Systems Mammomat Inspiration
Toshiba	Ultrasound Equipment
Tough-PAC	iPad Covers
Universal Ferrule Remover	Ferrule Remover
Zoll	AED

\*Previously Clinell Sporicidal Wipes.

Approved Manufacturers. © GAMA Healthcare Ltd. Version 14.

www.gamahealthcare.com/compatibility

**GAMA Healthcare Ltd.**, The Maylands Building, Maylands Avenue, Hemel Hempstead, Hertfordshire, HP2 7TG, UK. +44 (0)20 7993 0030 | info@gamahealthcare.com | www.gamahealthcare.com

